



Facial Massagers - Fyola

Facial Massager is proven effective natural skin care remedies. Skin treated with our revolutionary facial massager looks younger, because it is healthier from the inside out. This facial massager is similar to the expensive ultrasound machines used in beauty salons. Now you can have this beauty skin care product at your home!



Zoom In

Revolutionary Anti-Aging Product

Fyola Facial Massager - Benefits

- ✓ Helps Remove Furrow
- ✓ Prevent Acne
- ✓ Remove old cuticle, oiled waste out of pores
- ✓ Enhance absorption of moisture & nutrition
- ✓ Ultrasonic Waves + Galvanic Ion
- ✓ Light Weight & Portable (cordless)

Add to Cart

Information

Fyola Facial Massager

This advanced skin care device includes an artificial intelligence controller to help you remove skin waste by using ultrasonic waves. In addition, Fyola facial massager enhances the absorption nutrition deep into you skin by using galvanic ions.

Fyola has proven to create and maintain a clean and healthy skin. Fyola can also be used to control body fat by massaging the skin tissues.

Fyola facial massager has 6 different modes:

1. **Cleansing**
2. **Massaging**
3. **Nutrition**
4. **Lifting**
5. **Whitening**
6. **Wrinkle Care**

Specifications

Ultrasonic Frequency: 1 MHz

Galvanic Ion: YES

Rechargeable: YES

Wireless: YES

Treatment Area: 4 cm²

Power: AC 100V/220V, 50/60Hz

Size: 160(L) x 30(W) x 40 (D) mm

Weight: 200g



A. Your skin blocks the absorption of nourishment products



B. Fortunately, Fyola enhances the absorption of nourishment products



6 Functions Explained

Step 1 - Cleansing Mode

The cleansing mode of the facial massager removes skin waste, keratin or sebum. The cleansing mode of the Fyola Facial Massager gives you a clean and clear skin through Ionzyme by using galvanic ion and ultrasonic waves together.



The Cleansing mode forces keratin, sebum and other skin waste to be scattered up on the epidermis and dermis off the epidermal layer by dissolving them.

Step 2 - Massaging Mode

Massaging mode helps make skin tissues elastic and more beautiful. It stimulates the metabolism by using positive/negative ions and ultrasonic waves together. The massaging mode activates your skin cells.



Massaging your skin will make the skin more



A basic understanding of human skin

Human skin tissues consist of an outermost epidermis (epidermal layer), fibrous dermis (dermal layer) and adipose hypodermis (hypodermis layer) that regenerates every 28-day. The skin barrier is very important in order to protect your skin; however it also has a negative affect in respect to beauty treatments, because it interferes with the absorption of nourishment products such as creams and lotions.



Add to Cart

How Fyola (Ionzyme) can help your Skin?

Galvanic-ionizing

Fyola enhances the absorption of nutriments by bringing the nourishment products deep into your skin using the anode (+) and the cathode (-) as of electric currents. This will reduce the pimples on your face and enhance the skin whitening.

Ultrasonic Function

Fyola is a unique beauty skin care product, which activates your


healthy skin cells and tissues through a million of micro vibrations per second. Ultrasonic waves give you an elastic and healthy skin through its washing-away thermo-vibrations.

What are the benefits of Fyola (Ionzyme)?

"Ionzyme" is a new skin-care method, which removes skin wastes and enhances the absorption of nourishment products by using the electrophoresis method and ultrasonic waves.

By using "Ionzyme", Fyola facial massager can reduce pimples, freckles, skin dryness/oiliness and tighten drooped eye pouches in much shorter time than conventional facial massagers. It also contract hair pores so you can have a healthy and clean skin, along with a beautiful and elastic body figure.

Ionzyme involves the use of an electrical current to facilitate greater penetration of vitamins into the skin. With conventional iontophoresis, changes are very much slower, whereas with the Fyola Facial massager, you will notice that your skin looks and feels firm and tight after two or three treatments.



Why use a Gel with Facial Massager?

A gel is needed to conduct the ultrasonic waves deep into the skin. Research shows that the transmission proportion of ultrasonic waves in the air is 0%, about 50% in general water, 60% in distilled water, near 70% in glue matter, and above 70% in special conductive gel. Besides, another positive effect of ultrasonic wave is the promotion of skin absorbing nutrition. A gel enhances the transmission of ultrasonic waves into your skin.

For best results, apply conducting gel to the area to be treated (they must be certified by hygiene organization). Do NOT apply cream when using ultrasonic waves.

How to use the Facial Massager?

1. Apply gel or lotion

Apply a suitable amount of water soluble gel or lotion (a circle with 2 cm diameter) on the facial head or your skin.



elastic and healthy by activating skin cells and tissues on the epidermal and dermal layer through a quick stir of positive ions and negative ions, and therefore promoted the regeneration of basal cells.

Step 3 - Nutrition Mode

Use the nutrition mode to moisturize your skin by bringing the nutrition products deep into your skin tissues through ionic inducements



Step 4 - Lifting Tight

Lifting mode of the facial massager makes your skin tissues elastic and more beautiful. It activates the elastin and collagen by using positive/ negative ions and ultrasonic waves.



Use the lifting mode to reduce wrinkles and make your skin tissues more elastic. By using the lifting mode of the facial massager, it will tighten your sagging facial muscle by using positive and negative ions.

Step 5 - Whitening Mode

Whitening mode makes your skin more clear and sheer by bringing the nutrition products deep into the skin tissues through galvanic-ionic and ultrasonic waves.

Whitening evens skin tone, smoothes skin texture, reduces ruddiness, fades freckles.

Do not use cream!



2. Switch the power ON

3. Choose the right mode

- Cleansing
- Massaging
- Nutrition
- Lifting
- Whitening
- Wrinkle Care



Step 6 - Wrinkle Care Mode

Wrinkle care mode reduces the wrinkles of your face. It makes your skin cells tissues more elastic and beautiful by activating skin cells through ultrasonic vibrations and positive/negative ions (galvanic ion).

4 Select Ultrasonic or Galvanic Ion

Select ultrasonic wave mode to remove waste inside your skin or galvanic ion mode to enhance the absorption of ingredients (Vitamin C) into your skin.



When ultrasonic and galvanic ion mode are used at the same time, the permission through the out skin layer will be 77 times better than using the ultrasonic mode or galvanic mode separately.

3. Move the facial head around your face

Move the treatment head of the facial massager in circle following the formation of the muscles. See illustrations for the directions. Do not stop moving the treatment head in use, about the rate of 5 centimeters within 1 second. The power would shut off automatically.

4. Wash your face

After using the facial massager, clean your skin with water or tissue papers

6. Clean the Facial Massager

Clean the treatment head with tissue papers after use.

Safety & Cautions