

## Why gives Fyola Facial Massage


**Fyola Facial Massager is the world first wireless rechargeable face massager with the combined functions of ions and ultraviolet that has been published in the academic journals in Japan, the UK and the US, and has received an International Patent!**

### Not All Galvanic Ion / Ultrasonic Facial Massager Are The Same!

Fyola facial massager has successfully completed the skin tests and clinical experiments laboratories at the Maysei University and Hiroshima University in Japan. Both universities have a high reputation for their achievements in genetic engineering. Fyola has proven to bring the effective usage of the galvanic ion into the skin!

Results of the tests and experiments have found that adopting galvanic ions instead of using hand massage, increases the absorption rate of pro-vitamin C 40.3 times more for the outer skin and 30.7 times more for the inner skin

Test Results Fyola Facial Massager by Hiroshima University in Japan





The combination of ions and ultraviolet increases the absorption 77.6 times for the outer skin more and 37.9 times more for the inner skin than only using hand massage.

These results scientifically demonstrate that Fyola skin care products are far more effective compared with other skin care products.

Using ions and ultraviolet simultaneously, Fyola Face Massager provides the most effective way to rejuvenate your skin and reduce your wrinkles

## A joint study involved by Maysei University & Lotts



- Fyola Face Massager completed the ion effect test organized by Hiroshima University, Japan
- Fyola Face Massager completed a test drawing galvanic ions organized by Maysei University, Japan
- Fyola Face Massager completed published in an academic journal, Fragrance in Japan
- Fyola Face Massager completed published in an academic journal, ICAVS-2 in the UK

- Fyola Face Massager completed received a silver prize at the 32nd Geneva International New Technology and Inventions Competition

**Ultrasonic Wave is invisible for human eyes**

The ultrasound skin care product emits a deep micro-vibration so rapid, that cannot be seen or felt at the surface of your skin. So don't think your ultrasonic massager isn't working!

To so see the transmission of the ultrasonic waves, put some water at the metal head of the face massager. You will see the vibration of water, caused by the ultrasonic waves.