



How to Reduce your Wrinkles within 28 days?

The ultrasonic vibration of the facial micro massage is transmitted deep into your skin, so that it promotes metabolism and restore your skin! In **only 10-15 minutes** daily, you'll see and experience dynamic changes in your face: your skin will appear **more youthful, smooth and healthy in as soon as 3 days.**

The only Skin Care Massager that combines ultrasonic waves, ion & infrared!

Ultrasonic / Ion Facial Massagers

Ultrasonic Facial Massagers: New Microdermabrasion!

Reduce wrinkles, freckles, stains and double chin with ultrasonic facial massager. Skin treated with our revolutionary facial massager looks younger, because it is healthier from the inside out. This facial massager is similar to the expensive ultrasound microdermabrasion machines used in beauty salons. Now you can have this skin care product at your home!

Dione Skin Massager – Benefits

- Effective acne and wrinkle treatment
- Reduce Wrinkle
- Moisturize Skin
- Skin Firming
- Revolutionary Anti-Aging treatment
- Deep Cleansing
- Simple and Easy to use

Information

and pure by promoting absorption of nourishment of cosmetics.

Save Money spent at Beauty Salons

The performance of this model is similar to the expensive ultrasound machines used at most of the modern beauty salons. Now you can have this beauty product at your home, without spending thousands of dollars at your beauty salons.

It is effective on acne and wrinkle treatment by massaging the muscle with High Frequency Ultrasonic waves. Cleanse, revitalize, and rejuvenate your aging skin with the latest ultrasonic technology!

Make Skin Tight & Firm: Ultrasonic & Ion Combined with ultrasonic and negative ion, the facial massager keeps your skin beautiful and elastic. The ultrasonic facial massager protects your skin's health from pimple, discoloration on the face, dry skin and oily skin. This facial massager will help improve the appearance of skin tone, texture, elasticity, enlarged pores, fine lines and wrinkles. It's effective on specially aged skin.

How to get Partial Slimming?

This portable skin care equipment, makes your body slimmer and healthier through ultrasonic, elimination of waste matter in pores, buildup of skin resilience and adipose combustion with multiple functions using a million-time ultrasonic micro-vibration per second, galvanic ion and vibration.

The Ultrasonic Face and Body Massager can be applied to most of the body parts including face, neck, arms, breasts,

belly, waist, hips, buttocks, and legs. This product will help you achieve healthy, younger looking, younger feeling skin as well as help the slimming of the face and body.

Why use a Gel with Facial Massager?

A gel is needed to conduct the ultrasonic waves deep into the skin. Research shows that the transmission proportion of ultrasonic waves in the air is 0%, about 50% in general water, 60% in distilled water, near 70% in glue matter, and above 70% in special conductive gel. Besides, another positive effect of ultrasonic wave is the promotion of skin absorbing nutrition. A gel enhances the transmission of ultrasonic waves into your skin.

For best results, apply conducting gel to the area to be treated (they must be certified by hygiene organization). Do NOT apply cream when using ultrasonic waves

Specifications

Power: AC 100V/220V, 50/60Hz

Ultrasonic Frequency: 1 MHz

Galvanic Ions: YES

Far Infrared: NO

Size: 170(L) x 57(W) x 57(D) mm

Weight: 200g

Gel: Not included



Safety & Cautions

1. To clean, gently wipe the ultrasonic massager head with cloth moistened in water or alcohol
2. Use Ultrasonic Massager on areas of muscles, tendons and ligaments only.
3. Avoid striking the ultrasonic massager head.
4. Do not apply the device near the eyes nor on areas where the skin is thin
5. Apply the Ultrasonic Massager by moving it lightly in a spiral motion. Do not leave the device on the same spot for over 3 seconds. Move the Ultrasonic Massager to another spot immediately when sensation of warmth is reached.

How to use the facial massager

1. Connect the facial massager to power source.

2. Apply suitable amount of gel (a circle with 2 cm diameter) on the treatment head
3. Switch the power ON
4. Move the treatment head of the facial massager in circle following the formation of the muscles. See illustrations for the directions. Do not stop moving the treatment head in use, about the rate of 5 centimeters within 1 second. The power would shut off automatically.
5. After using the facial massager, clean your skin with water or tissue papers.
6. Clean the treatment head with tissue papers after use.